



UBE2A/B rabbit pAb antibody

| Catalog No : | Source: | Concentration : | Mol.Wt. (kD): |
|--------------|---------|-----------------|---------------|
| A23104 | Rabbit | 1 mg/ml | 17 kD |

| | |
|------------------------------|--|
| Applications | WB,IHC,ELISA |
| Reactivity | Human |
| Dilution | WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. |
| Storage | -20°C/1 year |
| Specificity | UBE1L Polyclonal Antibody detects endogenous levels of UBE1L protein. |
| Source / Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen | The antiserum was produced against synthesized peptide derived from human UBE1L. AA range:963-1012 |
| Uniprot No | P41226 |
| Alternative names | UBA7; UBE1L; UBE2; Ubiquitin-like modifier-activating enzyme 7; Ubiquitin-activating enzyme 7; D8; Ubiquitin-activating enzyme E1 homolog |
| Form | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Clonality | Polyclonal |
| Isotype | |
| Conjugation | |
| Background | ubiquitin like modifier activating enzyme 7(UBA7) Homo sapiens The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three c |
| Other | Gene_name: UBA7 ; Protein_name: Ubiquitin-like modifier-activating enzyme 7; Expression: B-cell, Eye, Platelet, |

Product Images

Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.