



TAF II p100 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A22245	Rabbit	1 mg/ml	87 kD
Applications	WB,IHC,ELISA		
Reactivity	Human,Rat		
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	TAF I p48 Polyclonal Antibody detects endogenous levels of TAF I p48 protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human TAF1A. AA range:301-350		
Uniprot No	Q15573		
Alternative names	TAF1A; TATA box-binding protein-associated factor RNA polymerase I subunit A; RNA polymerase I-specific TBP-associated factor 48 kDa; TAFI48; TATA box-binding protein-associated factor 1A; TBP-associated factor 1A; Transcription factor SL1;		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	TATA-box binding protein associated factor, RNA polymerase I subunit A(TAF1A) Homo sapiens This gene encodes a subunit of the RNA polymerase I complex, Selectivity Factor I (SLI). The encoded protein is a TATA box-binding protein-associated factor th		
Other	Gene_name: TAF1A ; Protein_name: TATA box-binding protein-associated factor RNA polymerase I subunit A; Expression: Embryo,Skin,		
Product Images	<input type="checkbox"/>		

Application Key:



W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.