



CaMKIIβ/γ rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A11556	Rabbit	1 mg/ml	62 kD
Applications	WB,IF,ELISA		
Reactivity	Human,Mouse,Rat		
Dilution	WB: 1/500 - 1/2000. IF: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.		
Storage	-20°C/1 year		
Specificity	CaMKIIβ/γ Polyclonal Antibody detects endogenous levels of CaMKIIβ/γ protein.		
Source / Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
Immunogen	The antiserum was produced against synthesized peptide derived from human CaMK2 beta/gamma. AA range:509-558		
Uniprot No	Q13555/Q13554		
Alternative names	CAMK2G; CAMK; CAMK-II; CAMKG; Calcium/calmodulin-dependent protein kinase type II subunit gamma; CaM kinase II subunit gamma; CaMK-II subunit gamma; CAMK2B; CAM2; CAMK2; CAMKB; Calcium/calmodulin-dependent protein kinase type II subunit bet		
Form	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
Clonality	Polyclonal		
Isotype			
Conjugation			
Background	calcium/calmodulin dependent protein kinase II gamma(CAMK2G) Homo sapiens The product of this gene is one of the four subunits of an enzyme which belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein k		
Other	Gene_name: CAMK2B/CAMK2G ; Protein_name: Calcium/calmodulin-dependent protein kinase type II subunit beta/gamma; Expression: Brain,Epithelium,Eye,		
Product Images			

Application Key:



W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

Species Cross-Reactivity Key:

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All
Species Expected

Trademarks

All product names and trademarks are the property of their respective owners.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support:

To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).

To call, write, fax, or email us, please visit www.aabsci.com, contact information will be displayed.