



## ADH7 rabbit pAb antibody

Catalog No :	Source:	Concentration :	Mol.Wt. (kD):
A10352	Rabbit	1 mg/ml	40 kD
<b>Applications</b>	WB,ELISA		
<b>Reactivity</b>	Human,Monkey		
<b>Dilution</b>	WB: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.		
<b>Storage</b>	-20°C/1 year		
<b>Specificity</b>	ADH7 Polyclonal Antibody detects endogenous levels of ADH7 protein.		
<b>Source / Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.		
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human ADH7. AA range:211-260		
<b>Uniprot No</b>	P40394		
<b>Alternative names</b>	ADH7; Alcohol dehydrogenase class 4 mu/sigma chain; Alcohol dehydrogenase class IV mu/sigma chain; Gastric alcohol dehydrogenase; Retinol dehydrogenase		
<b>Form</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.		
<b>Clonality</b>	Polyclonal		
<b>Isotype</b>			
<b>Conjugation</b>			
<b>Background</b>	alcohol dehydrogenase 7 (class IV), mu or sigma polypeptide(ADH7) Homo sapiens This gene encodes class IV alcohol dehydrogenase 7 mu or sigma subunit, which is a member of the alcohol dehydrogenase family. Members of this family metabolize a wide		
<b>Other</b>	Gene_name: ADH7 ; Protein_name: Alcohol dehydrogenase class 4 mu/sigma chain; Expression: Esophagus,Foreskin,PCR rescued clones,Stoma		

### Product Images

#### Application Key:

W-Western IP-Immunoprecipitation IHC-Immunohistochemistry ChIP-Chromatin Immunoprecipitation  
IF-Immunofluorescence F-Flow Cytometry E-P-ELISA-Peptide

**Species Cross-Reactivity Key:**

H-Human M-Mouse R-Rat Hm-Hamster Mk-Monkey Vir-Virus Mi-Mink C-Chicken Dm-D. melanogaster  
X-Xenopus Z-Zebrafish B-Bovine Dg-Dog Pg-Pig Sc-S. cerevisiae Ce-C. elegans Hr-Horse All-All  
Species Expected

**Trademarks**

*All product names and trademarks are the property of their respective owners.*

**Regulatory Disclaimer**

*For life science research only. Not for use in diagnostic procedures.*

---

**Contact and Support:**

*To ask questions, solve problems, suggest enhancements and report new applications, please visit our [Online Technical Support Site](#).*

*To call, write, fax, or email us, please visit [www.aabsci.com](http://www.aabsci.com), contact information will be displayed.*